

Elite™

Digital Visual Acuity Chart




**Any screen
is your
visual acuity
chart !**

- Fits any LCD monitor or TV set.
- No computer needed.
- Work out of the box, plug and play.
- Full-functional and multiple charts.
- RF wireless remote, fluorescent keypad.
- High definition, up to 1920 X 1080 (HD 1080P).



Elite™ Digital Visual Acuity Chart



- Fluorescent keypad, testing in dark room is much easier.
- No aiming! 2.4G RF remote control is omnidirectional.
- Longer working distance, faster response, reliable performance.
- Optional power supply solution, USB port or charger.
- HDMI or VGA port optional.
- Compact size.
- Video, audio and pictures, customized storage and display.
- Smart calibration, easy setup.
- Continuously testing distance, from 1m to 7m.
- Random sequencing of optotypes.
- Imperial and Metric standard optional.
- Mirror function available.
- Online upgrade.
- Package size : 180mm X 140 mm X 70mm

